

→  
REF. M1-ETA1  
EMERGENCY TROLLEY



### TECHNICAL FEATURES



Width:	67 cm
Depth:	50 cm
Height:	100 cm (surface)
Ø of the wheels:	10 cm
Weight:	52 kg
Max. Load:	50 kg / drawer
Material:	Aluminum alloy & ABS
Surface:	White ABS
Color:	Red

### SPECIFICATIONS



- Soft closing sliding rails
- Lockable wheels
- Central lock with key

### INCLUDES



- Defibrillator shelf
- IV holder
- Oxygen tank holder
- CPR board
- 5 drawers
- 2 waste bins
- Needle box
- 304 Stainless steel wire basket
- Socket base
- Extendable side tray
- Dividers to organize the drawer

