

→  
REF. EL35B  
ELECTRIC BLOOD  
EXTRACTION CHAIR



## TECHNICAL FEATURES



Length:	178 cm
Width:	59 cm
Height:	55,4-86,5 cm
Weight:	80 kg
Max. Load:	225 kg
Motors:	1
Foam thickness:	5 cm
Density:	36 kg/m <sup>3</sup>
Upholstery:	MED PVC ●

## MOTOR FUNCTIONS



- Height (flat position): 55,4-86,5 cm

## MANUAL FUNCTIONS



- Backrest (0° to +65°) inclination by gas lift
- Footrest (0° to -90°) inclination by gas lift
- Adjustable armrest

## INCLUDES



- Medical control box
- Foot controller
- 4 retractable wheels system for easier movement
- Adjustable foot (to stabilize)

## CERTIFICATIONS CE




- 93/42/EEC Medical Device Directive

# Therp

→  
TECHNICAL DATA SHEET

CONFY CARE

→  
REF. EL35B  
ELECTRIC BLOOD  
EXTRACTION CHAIR



**KG** 225 kg  
max

